


## Stainless Steel Carboy Shelves

## Pegboard mount carboy shelf



<i>Features</i>	<i>Model Number</i>	<i>Dimensions</i>
• All stainless steel construction	CS-18-P	18"w x 12"d
• Made from 18 gauge material	CS-20-P	20"w x 12"d
• Type 304 stainless steel, #4 finish	CS-24-P	24"w x 12"d
• 4"h wraparound	CS-30-P	30"w x 12"d
	CS-32-P	32"w x 12"d
	CS-36-P	36"w x 12"d
	CS-48-P	48"w x 12"d