


# MOD-RACK™ stainless steel pegboard

Model # V3024




## "V" Victoria Series

V series pegboards have a catch drain and integral drip trough in 2" or 4" depth. All drip troughs come with a 36" clear flexible drain tube.


## Features

- 30" w x 24" h x 3/4" d
- 25- 6" long 1/2" diameter white polypropylene pegs
- Autoclavable
- Constructed from 20 gauge, type 304 stainless steel, polished to a #4 finish.
- Wall mount bracket shipped with each pegboard.

## Standard Peg



## Product Dimensions



Wall Mount Bracket

