


FINISHED FUME HOOD DECK SIZES				
A	B	C	D	E
36"	30"	4.81"	3.50"	3.44"
47.75"	31"	3.96"	6.00"	3.71"
48"	30"	5.58"	2.70"	2.50"
48"	30"	5.58"	3.20"	2.00"
48"	30"	5.21"	3.00"	3.45"
48"	31"	5.58"	2.50"	3.70"
59.75"	31"	3.95"	6.00"	3.71"
60"	30"	6.06"	2.70"	2.50"
60"	30"	6.06"	3.20"	2.00"
60"	30"	5.21"	3.00"	3.45"
60"	31"	6.06"	2.50"	3.70"
71.75"	31"	3.93"	6.00"	3.71"
72"	30"	6.06"	2.70"	2.50"
72"	30"	6.06"	3.20"	2.00"
72"	30"	5.18"	3.00"	3.45"
72"	31"	6.06"	2.50"	3.70"
95.75"	31"	3.91"	6.00"	3.71"
96"	30"	6.03"	2.70"	2.50"
96"	30"	6.03"	3.20"	2.00"
96"	30"	5.16"	3.00"	3.45"
96"	31"	6.03"	2.50"	3.70"
48"	36"	4.88"	3.50"	3.62"
60"	36"	4.88"	3.50"	3.62"
72"	36"	4.88"	3.50"	3.62"
96"	36"	4.79"	3.50"	3.56"
96"	42"	4.88"	3.50"	3.63"

TABLE REPRESENTS STANDARD SIZES, ADDITIONAL SIZES ARE AVAILABLE UPON REQUEST.



<i>PII</i> PRIME industries, inc.		
www.piiab.com		
FUME HOOD DECK SIZES		
DATE:	10.18.13	REVISED: 10.18.13
		DRAWING NUMBER: FUME HOOD DECKS